SOLUTIONS INNOVANTES EN MEDECINE

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**High-flow Nasal
Cannula Oxygen
Therapy (HFNO, HFNC)**

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The high-oxygen therapy device HF60L is designed to provide respiratory therapy to patients with respiratory failure including viral pneumonia.

The HFNO high-throughput device nasal therapy method is recommended by the World Health Organization for the treatment of oxygen deficiency and acute respiratory disease COVID-19 caused by the coronavirus SARS-CoV-2.

The HF60L device delivers a mixture of oxygen, air and water vapor through a nasal cannula at a rate of 6 to 60 liters per minute.

The flow rate, humidity and heating level are electronically controlled. The level of oxygen supplied to the mixture is adjusted by the doctor based on the patient's oxygen saturation data, taking into account the pulse oximetry, arterial blood gases and acid-base metabolism parameters.

Steam is converted using high frequency heating technology: the frequency generated by the electronic microprocessor system is converted into mechanical vibration by a piezoceramic crystal.

Mechanical vibrations transform water into vapor, while the water is not heated.

The steam temperature is electronically controlled to maintain a comfortable temperature of the oxygen-air-water mixture.

The HF60L can be used simultaneously in high-flow oxygen therapy and in nebulizer mode, allowing drugs to be added to the mixture.

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Built in blower reaches a flow rate of up to 60 liters per minute

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Designed for continuous operation

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